

## **ecl-machine-loc**<sup>11,40</sup>

$\forall i:\text{Id}, ds:\text{fpf}(\text{Id}; x.\text{Type}), da:\text{fpf}(\text{Knd}; k.\text{Type}), A:\text{ecl}(ds; da), snd:\text{msg-spec}(ds; da),$   
 $upd:\text{update-spec}(ds; da), j:\text{Id}.$   
 $\text{msg-spec-loc}(snd; i)$   
 $\Rightarrow \text{squal}(\text{R-has-loc}(\text{ecl-machine}\{\text{ecl:ut2}\}(i; ds; da; A; snd; upd); j); \text{eq\_id}(i; j))$