

ecl-machine3-loc^{11,40}

```

∀ds:fpf(Id; x.Type), da:fpf(Knd; k.Type), x,T:top, ks:(top List), a:top, snd:msg-spec(ds;
                                                                    da),
  j:Id.
sqequal(R-has-loc(ecl-machine3(ds; da; x; T; ks; a; snd); j);
  if null(ks)
  then reduce((λl,b. bor(band((¬b.null(ecl-tags(l; snd))); eq_id(source(l); j)); b));
              ff;
              remove-repeats(idlnk-deq; msg-spec-links(snd)))
  else reduce((λl,b. bor(eq_id(source(l); j); b));
              ff;
              remove-repeats(idlnk-deq; msg-spec-links(snd)))
  fi )
```